

Rapigen

One Step SARS-CoV-2 Antigen Rapid Test

BIOCREDIT COVID-19 Ag Home Test Nasal

One step, rapid, immuno-chromatographic assay for the qualitative detection of SARS-CoV-2 antigen from Nasal swab

BIOCREDIT COVID -19 Ag Home Test Nasal is a rapid chromatographic immunoassay intended for the qualitative detection of nucleocapsid (N) protein antigen from SARS-CoV-2 in individuals with or without symptoms or other epidemiological reasons to suspect COVID-19 infection.

This test is authorized for home use with **self-collected** direct nasal (nares) samples using the sterilized nasal swab provided in the kit.

1st and officially approved by MFDS

- 1 Quick test result within 15-30 mins
- 2 Easy to use on site
- 3 Accurate test
- 4 User friendly Dual color system
- 5 Mutant SARS-CoV-2 detection



MFDS **IVD**

No compromise on quality

BIOCREDIT COVID-19 Ag Home Test Nasal


One step, rapid, immuno-chromatographic assay for the qualitative detection of SARS-CoV-2 antigen from Nasal swab

Assay Procedure


1 Preparation before test




2 Assay diluent solution




3 Swab




4 Filter cap



5 Test device




Interpretation of Results



Description of materials provided

Sterilized nasal swab, Filter cap, Assay buffer, Test device, Instructions for use, Plastic bag for disposal.



Manufactured by —

Rapigen, INC.

1152SNH073A00

161, Saneop-ro 155beon-gil, Gwonseon-gu, Suwon-si, Gyeonggi-do 16648, Republic of Korea
TEL : +82-31-427-4677 FAX : +82-31-427-4678 E-mail : info@rapigen-inc.com www.rapigen-inc.com